

Chromacity 520

520 nm ultrafast source



01 Applications


- Pump source for OPO
- Non-linear optics
- Telecoms / quantum imaging
- Raman spectroscopy
- Photo-dynamic therapy
- Nanophotonics

02 Technical Overview

- Pulse duration < 150 fs
- 1.5 W average power available
- Access to 1040 nm pulses

03 Features & Benefits

- Compact laser housing
- Intuitive web browser interface
- Power efficient laser cavity (low power consumption)







Pumped by a Chromacity 1040 ensures delivery of high quality ultrashort pulses in the green.

Chromacity 520

520 nm ultrafast source

Specifications	Low Power Chromacity 1040 and internal SHG module	High Power Chromacity 1040 and external SHG module
Output power	Up to 500 mW	Up to 1.5 W
Wavelength	520 nm (1040 nm also available)	
Pulse duration	< 150 fs	
Repetition frequency	100 MHz	
Beam parameters	Free space, Linearly polarized Divergence < 2 mrad	
Beam diameter	0.5-1.2 mm	
Control interface	Web browser interface Ethernet & serial port (RS232) also available	
Electrical	Voltage 110 – 240 V AC, Frequency 50 – 60 Hz, Power 80 W	
Dimensions	PSU 483 x 255 x 177 mm	
	Laser head 526 x 250 x 108 mm	Laser head 426 x 250 x 108 mm SHG module 150 x 300 mm



Chromacity follows a policy of continuous improvement, hence specifications are subject to change without notice.

Learn how our ultrafast lasers can enable you to discover more.
For more information, email: sales@chromacitylasers.com



Discover More
chromacitylasers.com